




LCMS Report

Sample Name : M5680-Gemcitabine-03
Method File : B 5%-95%-P.lcm

Method info. : Waters SunFire C18 50*4.6mm 5um 2.000ml/min 2.6min Column Temperature: 40 °C
Gradient: 5% B hold for 0.2min, increase to 95 % B within 1.40 min, hold at 95 % B for 0.9 min, then back to 5% B within 0.01 min

Pump A : 0.1% FA in H2O
Pump B : 0.1% FA in ACN



Integration Results for DAD

PDA Ch1 254nm

Peak#	Ret. Time	Height	Area	Area%
1	0.400	782324	1692419	100.000
Total		782324	1692419	100.000

PDA Ch2 214nm

Peak#	Ret. Time	Height	Area	Area%
1	0.402	1273132	2785197	100.000
Total		1273132	2785197	100.000

Mass Spectrum

Rt=0.403min Positive

