




Sample name: M5785-Metoprolol Tartrate-01  
Instrument: 1.5  
Inj. volume: HPLC  
Acq. method:  
Processing method:  
Manually modified:

Injection date:  
Location:  
Type:  
Sample amount:



Signal: DAD1A,Sig=220,4 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
1.890		34.0	11.2	0.901
5.355		0.6	0.2	0.017
5.574		3741.5	1021.9	99.082
	Sum	3776.21		