



Sample name: M5929-Prostaglandin E2-05

Injection date:

Instrument: 3

Location:

Inj. volume: HPLC


Type:

Acq. method:

Sample amount:

Processing method:

Manually modified:



Signal: DAD1C,Sig=254,4 Ref=off

RT [min]	Peak Resolution USP	Width [min]	Area	Height	Area%
5.870		0.05	35.2	13.0	0.22
6.533		0.16	5.5	1.0	0.03
6.983		0.07	15541.1	3532.1	98.09
7.138		0.03	92.1	46.7	0.58
7.286		0.03	61.1	31.4	0.39
7.424		0.04	7.0	3.1	0.04
7.590		0.04	70.3	30.8	0.44
7.765		0.06	32.0	9.6	0.20
		Sum	15844.23		