




Sample name: M6078-CFI-402257-03  
Instrument: 2  
Inj. volume: HPLC  
Acq. method:  
Processing method:  
Manually modified:

Injection date:  
Location:  
Type:  
Sample amount:



Signal: DAD1D, Sig=235,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%	Name
5.424	MM m	0.04	1.13	0.41	0.01	
5.881	MM m	0.06	1.26	0.31	0.02	
6.044	MM m	0.04	0.27	0.10	0.00	
6.097	MM m	0.04	0.19	0.08	0.00	
6.294	BB	0.42	8147.93	2169.70	99.90	
6.626	MM m	0.04	0.79	0.29	0.01	
6.859	MM m	0.05	1.93	0.58	0.02	
6.971	MM m	0.05	0.77	0.26	0.01	
7.097	MM m	0.05	1.54	0.51	0.02	
	Sum		8155.79			