




Sample name: M6118-666-15-03
Instrument: 1.5
Inj. volume: HPLC
Acq. method:
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1H,Sig=240,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%	Name
6.708	MM m	0.05	1.36	0.47	0.01	
7.083	MB m	0.06	4.88	1.33	0.03	
7.656	MM m	0.05	0.93	0.29	0.01	
7.878	MM m	0.05	72.03	22.59	0.47	
8.483	BV	0.59	15222.55	3329.34	99.31	
9.495	MB m	0.06	26.56	6.36	0.17	
	Sum		15328.32			