




Sample name: M6199-MSC2530818-01  
Instrument: A010  
Inj. volume: 1.5  
Acq. method: HPLC  
Processing method:  
Manually modified:

Injection date:  
Location:  
Type:  
Sample amount:



Signal: DAD1D,Sig=230,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
5.942	MV m	0.13	4.39	0.97	0.04
6.944	BB	0.19	48.79	15.07	0.41
7.644	BB	0.63	11742.91	2942.67	99.16
8.200	BB	0.18	20.92	5.93	0.18
8.612	MB m	0.17	17.25	4.56	0.15
9.661	MB m	0.16	8.31	1.95	0.07