




Sample name: M6286-Selonserib-01  
Instrument: 2  
Inj. volume: HPLC  
Acq. method:  
Processing method:  
Manually modified:

Injection date:  
Location:  
Type:  
Sample amount:



Signal: DAD1D,Sig=230,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
5.377	MM m	0.07	0.42	0.19	0.00
5.448	MM m	0.09	0.76	0.29	0.01
5.590	BB	0.11	6.41	2.29	0.06
5.685	BB	0.12	5.88	2.09	0.05
5.903	BM m	0.22	10703.64	2939.73	98.74
6.093	MM m	0.13	89.84	27.50	0.83
6.240	BB	0.17	27.13	8.80	0.25
7.123	MB m	0.11	5.62	1.78	0.05