




Sample name: M8400-Fialuridine-02
Instrument: A010
Inj. volume: 1.5
Acq. method: HPLC
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1A, Sig=220.0, 4.0 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
4.507		6278.6	1946.1	99.817
6.866		5.2	1.6	0.083
11.904		2.5	0.6	0.040
12.220		2.0	0.5	0.032
12.344		1.7	0.4	0.028