




Sample name: M8992-AM966-01
Instrument: A010
Inj. volume: 1.5
Acq. method: HPLC
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1B, Sig=210,4 Ref=off

RT [min]	Area	Height	Peak	Area%
7.555	3.83	1.21		0.05
7.954	4.57	1.50		0.06
8.881	6.48	1.86		0.08
8.960	2.74	0.83		0.04
9.259	12.77	4.08		0.17
10.013	7668.70	1328.89		99.61