




Sample Name : M9243-ITD-1-01
Acq. Operator : C22
Spec. Reported : MS Integration
Acq. Method : D:\METHODS\2-POS-MON-1.M
Analysis Method : D:\METHODS\2-POS-MON-1.M
Sample Info : Easy-Access Method: '2-POS-MON-1'
Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)
Gradient: 5%-95% B in 1.5min, 95%B for 1.5 min
Flow Rate: 2ml/min
Column: Sunfire, 4.6*50mm, 3.5um
Oven Temperature: 50 C
MASS Range: 120-1000



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
0.51	0.01	45.60	64.11	1.44
0.54	0.02	10.22	9.60	0.32
2.28	0.02	3110.77	2460.87	98.24

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
2.28	0.02	3317.65	2797.35	100.00

Ret. Time: 2.30 <<<< POSITIVE SPECTRA >>>>

