




Sample name: M9428-MRTX849-04
Instrument: 1
Inj. volume: HPLC
Acq. method:
Processing method:
Manually modified:

Injection date:
Location:
Type:
Sample amount:



Signal: DAD1A,Sig=220,4 Ref=off

RT [min]	Peak Resolution USP	Area	Height	Area%
5.728		1.7	0.5	0.011
6.295		2.9	1.0	0.019
6.706		1.9	0.5	0.012
6.850		27.7	5.4	0.183
7.249		15078.1	2871.0	99.345
7.553		16.3	4.4	0.108
7.718		13.6	4.0	0.090
8.054		13.2	3.0	0.087
8.200		14.0	3.1	0.093
8.532		8.1	1.9	0.053
	Sum	15177.60		