




Agilent Technologies

Sample name: M9477-NSC348884-01  
Instrument: 1  
Inj. volume: HPLC  
Acq. method:  
Processing method:  
Manually modified:

Injection date:  
Location:  
Type:  
Sample amount:



Signal: DAD1B,Sig=210,4 Ref=off

RT [min]	Area	Height	Area%
5.535	14.86	3.74	0.09
6.256	11.61	1.77	0.07
6.683	12.30	4.05	0.07
6.857	16833.06	1225.41	99.76
8.013	2.12	0.89	0.01
Sum	16873.95		