




LC/MS Report

Injection Date : 31 Mar 21 :
 Sample Name : SP-13786-03
 Method Info : Column:SunFire C18 (4.6x 50 mm,3.5 um)
 Mobile phase: H2O(0.01%TFA) (A) / ACN(0.01%TFA) (B)
 Elution program: Gradient from 10 to 95% of B in 1.4 min at 2ml/min
 Temperature: 50 °C
 Detection: UV (214 , 4 nm) and MS (ESI, Pos mode ,110 to 1000 amu)



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.41	0.02	700.50	399.59	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.41	0.03	81.93	46.39	100.00

Ret. Time: 1.41

<<<< POSITIVE SPECTRA >>>>

