



Sample name: M9713-SJ000291942-01

Injection date:

Instrument: A010

Location:

Inj. volume: 3


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1D,Sig=230.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
8.037	MB m	0.16	12.47	3.79	0.06
9.827	MB m	0.20	8.98	2.04	0.04
10.960	BV	0.65	20160.67	2781.86	99.63
11.160	VB	0.11	7.79	1.81	0.04
11.674	MB m	0.19	34.30	8.07	0.17
12.757	MB m	0.16	10.47	2.45	0.05