




Sample Name : M9794-Ebselen-01
Acq. Operator : C22
Spec. Reported : MS Integration
Acq. Method : D:\METHODS\2-POS-MON-1.M
Analysis Method : D:\METHODS\2-POS-MON-1.M
Sample Info : Easy-Access Method: '2-POS-MON-1'
Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)
Gradient: 5%-95% B in 1.5min, 95%B for 1.5 min
Flow Rate: 2ml/min
Column: Sunfire, 4.6*50mm, 3.5um
Oven Temperature: 50 C
MASS Range: 120-1000



Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
1.19	0.01	21.23	22.82	0.87
1.26	0.01	20.64	27.01	0.84
1.93	0.02	2389.69	2030.48	97.78
2.11	0.02	12.29	9.91	0.50

Integration Results for DAD1 B, Sig=254,4 Ref=off


RetTim	Width	Area	Height	Area%
1.93	0.02	1156.88	1013.62	99.12
2.11	0.02	10.30	8.07	0.88

Ret. Time: 0.38

<<<< POSITIVE SPECTRA >>>>



Ret. Time: 1.28



Ret. Time: 1.95

