


Current Data Parameters
 NAME ONX 0914
 EXPNO 1
 PROCNO 1

F2 - Acquisition Parameters
 Date_ 20140530
 Time 18.52
 INSTRUM spect
 PROBHD 5 mm PABBO BB-
 PULPROG zg30
 TD 65536
 SOLVENT DMSO
 NS 16
 DS 2
 SWH 8278.146 Hz
 FIDRES 0.126314 Hz
 AQ 3.9584243 sec
 RG 362
 DW 60.400 usec
 DE 8.00 usec
 TE 297.3 K
 D1 1.00000000 sec
 TD0 1

===== CHANNEL f1 =====
 NUC1 1H
 P1 10.70 usec
 PL1 -3.00 dB
 SFO1 400.1324710 MHz

F2 - Processing parameters
 SI 32768
 SF 400.1300034 MHz
 WDW EM
 SSB 0
 LB 0.30 Hz
 GB 0
 PC 1.00



1.006
 1.000
 0.977
 5.543
 2.155
 2.160
 1.077
 1.068
 1.092
 3.287
 4.274
 0.667
 1.048
 1.060
 0.486
 2.776
 0.760
 0.348
 1.248
 0.988
 4.564
 3.199
 3.425

9 8 7 6 5 4 3 2 1 0 ppm