


LCMS Report

mAU




Peak Table

PDA Ch1 220nm

Peak#	Ret. Time	Height	Height%	Area	Area%
1	1.187	3454573	100.000	14389641	100.000
Total		3454573	100.000	14389641	100.000


mV



Peak Table


ELSD

Peak#	Ret. Time	Height	Height%	Area	Area%
1	1.184	1023239	100.000	3878515	100.000
Total		1023239	100.000	3878515	100.000



Peak Table

Peak#	Ret. Time	Height	Height%	Area	Area%
1	1.179	1747228	100.000	4854442	100.000
Total		1747228	100.000	4854442	100.000



Peak Table


Peak#	Ret. Time	Height	Height%	Area	Area%
1	1.184	1023239	100.000	3878515	100.000
Total		1023239	100.000	3878515	100.000

Segment#1 (x1,000,000)

2,188,233
(-0.055min Delay)

2.00

1.00



Retention Time: 1.178

Spectrum Mode:Single 1.178(149) Base Peak:441.30(533953)

BG Mode:None Segment 1 - Event 1

