


Injection Date : 29 Oct 20
Sample Name SAR407899-02
Method Info : Mobile Phase : A:water(0.01%TFA) B:ACN(0.01%TFA)
Gradient:5% increase to 95%B within 1.3min,95%B for 1.7min
Flow Rate :2ml/min
Column :Sunfire, 50*4.6mm,3.5um
Column Temperature:50C
Detection:UV(214,4nm) and MS(ESI,Pos mode,110 to 1000 amu)



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
--------	-------	------	--------	-------

Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.15	0.02	967.60	673.86	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.15	0.02	1951.81	1340.00	100.00

Integration Results for MSD1 TIC, MS File

RetTim	Width	Area	Height	Area%
1.15	0.04	1529438.88	593319.69	100.00

Ret. Time: 1.15 <<<< POSITIVE SPECTRA >>>>

