



Current Data Parameters
NAME SU11274
EXPNO 1
PROCNO 1

F2 - Acquisition Parameters

Date_ 20140531
Time 12.12
INSTRUM spect
PROBHD 5 mm PABBO BB-
PULPROG zg30
TD 65536
SOLVENT DMSO
NS 16
DS 2
SWH 8278.146 Hz
FIDRES 0.126314 Hz
AQ 3.9584243 sec
RG 256
DW 60.400 usec
DE 8.00 usec
TE 297.6 K
D1 1.00000000 sec
TD0 1

===== CHANNEL f1 =====

NUC1 1H
P1 10.70 usec
PL1 -3.00 dB
SFO1 400.1324710 MHz

F2 - Processing parameters

SI 32768
SF 400.1300000 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00

