


File D:\DATA\
Injection Date : 19 Nov 16 12:01 pm +0800 Tgt Mass (EZX) :
Sample Name : Saracatinib-01 Location : P1-C-02
Acq. Operator : LCMS Inj : 1
Spec. Reported : MS Integration Inj Volume : 2 ul
Acq. Method : D:\METHODS\1-POS-3MIN.M
Analysis Method : D:\METHODS\1-POS-3MIN.M
Sample Info : *** No target masses specified ***
Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)
Gradient: 5%-95% B in 1.2min
Flow Rate: 2.2 ml/min
Column: Poroshell 120 EC-C18, 4.6*30mm, 2.7um
Oven Temperature: 50 C
Mass Range: 110-1500



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
--------	-------	------	--------	-------


Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.30	0.07	1778.48	369.20	100.00


Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.30	0.06	1543.61	322.97	100.00


Ret. Time: 1.19



Ret. Time: 1.21



Ret. Time: 1.24



Ret. Time: 1.31

