



File D:

Injection Date : 21 Dec 16 07:07 pm
Sample Name : Sildenafil citrate-02
Acq. Operator : XXXX
Spec. Reported : MS Integration Inj Volume : 15 ul
Acq. Method : C:\Chem31\METHODS\1-POS-2MIN.M
Analysis Method : C:\Chem31\METHODS\1-POS-2MIN.M
Sample Info : Easy-Access Method: '1-POS-2MIN'
Method Info : Mobile Phase: A: Water(0.01%TFA) B: ACN(0.01%TFA)
Gradient: 5%B increase to 100%B within 0.8min,100%B for 1.1min.
Flow Rate:1.5mL/min
Column: Chromolith Fast Gradien RP-18e,50-3mm.
Column Temperature: 40C
Mass Range:100-1100


*ADC1 A, ELSD




*DAD1 A, Sig=214,4 Ref=off



*DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70,



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
0.85	0.03	211.66	105.89	100.00

Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
0.83	0.02	2119.31	1360.17	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
0.83	0.02	1026.62	667.78	100.00

Ret. Time: 0.86

<<<< POSITIVE SPECTRA >>>>

