

LCMS REPORT

Sample Information


LCMS

Acquired by : System Administrator Date Acquired : 2021/5/15 13:12:44
 Sample Name : Sorafenib Tosylate-07 injection Volume : 1
 Method File : B 5%-95%-P-new - 2.6min.lcm Vial : 10


Column : Waters SunFire C18 50*4.6mm 5um 2.000ml/min 2.6min
 Pump A : 0.1% FA and 10% ACN in H2O
 Pump B : 0.1% FA and 10% H2O in ACN

Chromatogram


mAU



mAU



段号1



Integration Results for DAD

PDA Ch1 254nm

Peak#	Ret. Time	Height	Area	Area%
1	1.895	3753355	6143201	100.000
总计		3753355	6143201	100.000


PDA Ch2 214nm

Peak#	Ret. Time	Height	Area	Area%
1	1.895	2914863	5220269	100.000


Peak#	Ret. Time	Height	Area	Area%
总计		2914863	5220269	100.000

Mass Spectrum


Rt=0.130min 负




Rt=0.259min 正




Rt=0.462min 正




Rt=0.805min 负




Rt=1.182min 正




Rt=1.840min 负




Rt=1.857min 正



Rt=2.155min 负



Rt=2.267min 负



Rt=2.262min 正

