



Current Data Parameters
NAME Tomeglovir-03
EXPNO 1
PROCNO 1

F2 - Acquisition Parameters
Date_ 20190424
Time_ 12.15
INSTRUM spect
PROBHD 5 mm PADUL 13C
PULPROG zg30
TD 65536
SOLVENT DMSO
NS 16
DS 2
SWH 8012.820 Hz
FIDRES 0.122266 Hz
AQ 4.0894465 sec
RG 199.46
DW 62.400 usec
DE 6.50 usec
TE 294.6 K
D1 1.00000000 sec
TD0 1

==== CHANNEL f1 =====
SFO1 400.1324710 MHz
NUC1 1H
P1 10.00 usec
PLW1 20.00000000 W

F2 - Processing parameters
SI 65536
SF 400.1300000 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00

Integration values (from left to right): 0.38, 0.98, 1.00, 2.03, 1.00, 0.42, 0.46, 2.48, 0.50, 2.04, 0.99, 2.07, 0.99, 2.01, 6.09, 6.06