

File D:\DATA\T55\MONITOR\

Injection Date : 29 May 17 12:09

Sample Name : Tubacin-01

Acq. Operator : LCMSA008

Spec. Reported : MS Integration

Acq. Method : C:\Chem32\1\METHODS\1-POS-3MIN.M

Analysis Method : C:\Chem32\1\METHODS\1-POS-3MIN.M

Sample Info : *** No target masses specified ***

Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)

Gradient: 5%-95% B in 1.7 min

Flow Rate: 2.2 ml/min

Column: SunFire C18, 4.6*50mm, 3.5um

Column Temperature: 50 C

MASS Range:110-900


Tgt Mass (EZX) :

Location : P1-F-02


Inj : 1

Inj Volume : 2 ul


*ADC1 A, ELSD




*DAD1 A, Sig=214,4 Ref=off



*DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Scan, Frag: 70



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
--------	-------	------	--------	-------

Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
2.03	0.02	895.84	664.41	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
2.03	0.02	722.86	559.37	100.00

LCMSA008 2017-05-29

Ret. Time: 0.30 <<<< POSITIVE SPECTRA >>>>



Ret. Time: 2.03

