



NAME vx950

EXPNO 1
PROCNO 1
Date_ 20110729
Time_ 14.44
INSTRUM spect
PROBHD 5 mm PABBO BB-
PULPROG zg30
TD 65536
SOLVENT CDCl3
NS 8
DS 2
SWH 8223.695 Hz
FIDRES 0.125483 Hz
AQ 3.9846387 sec
RG 203
DW 60.800 usec
DE 6.50 usec
TE 299.1 K
D1 1.00000000 sec
TD0 1

===== CHANNEL f1 =====

NUC1 1H
P1 13.00 usec
PL1 -2.80 dB
PL1W 18.85802460 W
SFO1 400.1324710 MHz
SI 32768
SF 400.1300097 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00

