


File D:\DATA\
Injection Date : 20 May 20 10:27 am +0800 Tgt Mass (EZX) :
Sample Name : AG-014699 phosphate-01 Location : P1-C-07
Acq. Operator : LCMS Inj : 1
Spec. Reported : MS Integration Inj Volume : 2 ul
Acq. Method : C:\Chem32\1\METHODS\1-POS-3MIN.M
Analysis Method : C:\Chem32\1\METHODS\1-POS-3MIN.M
Sample Info : Easy-Access Method: '1-POS-3min'
Method Info : Mobile Phase: A: water (0.01%TFA) B: ACN (0.01%TFA)
Gradient: 5%-95% B in 1.7 min
Flow Rate: 2.2 ml/min
Column: SunFire C18, 4.6*50mm, 3.5um
Column Temperature: 50 C
MASS Range: 110-900



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
--------	-------	------	--------	-------


Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.29	0.02	2125.34	1462.09	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.29	0.02	1178.35	831.66	100.00

Ret. Time: 1.31 <<<< POSITIVE SPECTRA >>>>



Ret. Time: 2.73

