

File D:\DATA\

Injection Date : 19 Nov 16 11:25 am +0800

Sample Name : AM2201-01

Acq. Operator : LCMS

Spec. Reported : MS Integration

Acq. Method : D:\METHODS\1-POS-3MIN.M

Analysis Method : D:\METHODS\1-POS-3MIN.M

Sample Info : *** No target masses specified ***

Method Info : Mobile Phase: A: water (10mM Ammonium hydrogen carbonate) B: ACN

Gradient: 5%-95% B in 1.5min

Flow Rate: 2.0ml/min

Column: XBridge C18, 4.6*50mm, 3.5um

Oven Temperature: 40 C

MASS Range: 110-800


Tgt Mass (EZX) :

Location : P2-A-04


Inj : 1

Inj Volume : 2 ul


*ELS1 A, Voltage




*DAD1 A, Sig=214,4 Ref=off



*DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70



Integration Results for ELS1 A, Voltage

RetTim	Width	Area	Height	Area%
1.51	0.05	363.58	119.89	100.00

Integration Results for DAD1 A, Sig=214,4 Ref=off



RetTim	Width	Area	Height	Area%
1.48	0.01	1372.90	1614.70	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.48	0.01	211.43	248.61	100.00

Ret. Time: 1.46 <<<< POSITIVE SPECTRA >>>>

Ret. Time: 2.29



Ret. Time: 2.40

Ret. Time: 2.48

