



File D:\DATA\
 Sample Name : AMG-458-01
 Acq. Method : D: HPLC.M
 Analysis Method : D: HPLC.M
 Sample Info : Easy-Access Method: '1-HPLC.M'
 Method Info : Mobile Phase: A: Water(0.01%TFA) B: ACN(0.01%TFA)
 Gradient: 5%B increase to 95%B within 8min,
 95%B for 7min, back to 5%B within 0.01min.
 Flow Rate: 1 mL/min
 Column: SunFire C18, 4.6*150mm, 3.5um
 Column Temperature: 40 C



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
6.15	0.04	2536.89	1000.58	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
6.15	0.04	2501.38	984.59	100.00