




LC/MS Report

Injection Date : 3 Mar 21
Sample Name : AMPPD-01
Method Info : Column:Xbridge C18(2) (4.6x 50 mm, 3.5um)
Mobile phase: H2O(10 mmol NH4HCO3) (A) / ACN(B)
Elution program: Gradient from 10 to 95% of B in1.5min at 1.8ml/min
Temperature: 50 °C
Detection: UV (214 , 4 nm) and MS (ESI, Pos mode ,103 to 800 amu)



Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
1.17	0.02	3076.25	1962.42	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off


RetTim	Width	Area	Height	Area%
1.17	0.02	163.67	108.55	100.00

Ret. Time: 1.16


<<<< POSITIVE SPECTRA >>>>



Ret. Time: 1.57



Ret. Time: 1.90



Ret. Time: 2.02

