

File D:\DATA

Injection Date : 28 Jul 16 06:45 pm

Sample Name : AN2728-01

Location :

Acq. Operator : XXXX

Inj : 1

Spec. Reported : MS Integration

Inj Volume : 10 ul

Acq. Method : D:\METHODS\1-POS-3MIN.M

Analysis Method : D:\METHODS\1-POS-3MIN.M

Sample Info : Easy-Access Method: '1-POS-3MIN'

Method Info : Mobile Phase: A: water(0.01%TFA) B:ACN(0.01%TFA)


Gradient: 5% to 95%B within 1.4 min

Flow Rate :2.0ml/min


Column :SunFire C18, 4.6*50mm,3.5um

Oven Temperature : 50C


*ADC1 A, ADC1 CHANNEL A




DAD1 A, Sig=214,4 Ref=off



DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70




Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.68	0.03	285.13	142.43	100.00


Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.69	0.03	276.00	136.28	100.00


Ret. Time: 0.85 <<<< POSITIVE SPECTRA >>>>




Ret. Time: 1.70



Ret. Time: 1.78



Ret. Time: 1.93



Ret. Time: 2.36

