




Injection Date : 26 Jan 21  
Sample Name : AZD3965-03  
Method Info : Column:SunFire C18 (4.6x 50 mm,3.5 um)  
Mobile phase: H2O(0.01%TFA) (A) / ACN(0.01%TFA)(B)  
Elution program: Gradient from 10 to95% of B in 1.4 min at 2ml/min  
Temperature: 50 °C  
Detection: UV (214 , 4 nm) and MS (ESI, Pos mode ,110 to 1000 amu)



Integration Results for DAD1 A, Sig=214,4 Ref=off

| RetTim | Width | Area    | Height  | Area%  |
|--------|-------|---------|---------|--------|
| 1.61   | 0.02  | 1341.39 | 1025.58 | 100.00 |


Integration Results for DAD1 B, Sig=254,4 Ref=off

| RetTim | Width | Area   | Height | Area%  |
|--------|-------|--------|--------|--------|
| 1.61   | 0.02  | 475.41 | 363.70 | 100.00 |

Ret. Time: 1.61

&lt;&lt;&lt;&lt; POSITIVE SPECTRA &gt;&gt;&gt;&gt;

\*ES-API Positive



Ret. Time: 1.77

\*ES-API Positive

