




11.870
11.840
11.827
11.791
11.753
11.695
11.673
11.620
8.257
8.253
8.249
8.241
8.235
8.223
8.198
8.192
8.184
8.179
8.176
8.176
7.966
7.962
7.948
7.943
7.941
7.933
7.925
7.923
7.918
7.904
7.899
7.691
7.670
7.251
7.231
3.342
3.328
3.235
2.514
2.510
2.506
1.235
0.014
0.006

Current Data Parameters
NAME Alizarin-01
EXPNO 1
PROCNO 1

F2 - Acquisition Parameters
Date_ 20190613
Time_ 12.37
INSTRUM spect
PROBHD 5 mm PADUL 13C
PULPROG zg30
TD 65536
SOLVENT DMSO
NS 16
DS 2
SWH 8012.820 Hz
FIDRES 0.122266 Hz
AQ 4.0894465 sec
RG 199.46
DW 62.400 usec
DE 6.50 usec
TE 297.7 K
D1 1.00000000 sec
TD0 1

===== CHANNEL f1 =====
SFO1 400.1324710 MHz
NUC1 1H
P1 10.00 usec
PLW1 20.00000000 W

F2 - Processing parameters
SI 65536
SF 400.1300000 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00



1.44
2.01
2.02
1.01
1.00