



File D:\DATA
Injection Date : 28 Jul 16 06:45 pm
Sample Name : Amuvatinib-02
Acq. Operator : XXXX
Spec. Reported : MS Integration
Acq. Method : D:\METHODS\1-POS-3MIN.M
Analysis Method : D:\METHODS\1-POS-3MIN.M
Sample Info : Easy-Access Method: '1-POS-3MIN'
Method Info : Mobile Phase: A: water(0.01%TFA) B:ACN(0.01%TFA)
Gradient: 5% to 95%B within 1.4 min
Flow Rate :2.0ml/min
Column :SunFire C18, 4.6*50mm,3.5um
Oven Temperature : 50C


*ADC1 A, ADC1 CHANNEL A




DAD1 A, Sig=214,4 Ref=off



DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70



Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
1.41	0.10	21180.91	2554.57	99.97
2.29	0.03	5.52	3.08	0.03

Integration Results for DAD1 B, Sig=254,4 Ref=off


RetTim	Width	Area	Height	Area%
1.41	0.11	24347.06	2957.33	100.00

Ret. Time: 0.32


<<<< POSITIVE SPECTRA >>>>




Ret. Time: 0.49




Ret. Time: 0.64




Ret. Time: 0.85




Ret. Time: 1.43




Ret. Time: 1.88




Ret. Time: 1.91




Ret. Time: 2.29




Ret. Time: 2.33




Ret. Time: 2.37




Ret. Time: 2.44




Ret. Time: 2.66



Ret. Time: 2.71



Ret. Time: 2.87



Ret. Time: 2.94

