




Current Data Parameters

NAME BAY-1895344-01
EXPNO 4
PROCNO 1

F2 - Acquisition Parameters

Date_ 20170703
Time 9.56
INSTRUM spect
PROBHD 5 mm PABBO BB/
PULPROG zg30
TD 65536
SOLVENT DMSO
NS 16
DS 2
SWH 8012.820 Hz
FIDRES 0.122266 Hz
AQ 4.0894465 sec
RG 131.1
DW 62.400 usec
DE 6.50 usec
TE 295.2 K
D1 1.00000000 sec
TD0 1

===== CHANNEL f1 =====

SFO1 400.1324710 MHz
NUC1 1H
P1 10.00 usec
PLW1 16.00000000 W

F2 - Processing parameters

SI 65536
SF 400.1300020 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00