



File D:\DATA\
Sample Name : CPI-613-04
Spec. Reported : MS Integration Inj Volume : 2 ul
Acq. Method : D: HPLC
Analysis Method : D: HPLC
Sample Info : Easy-Access Method:
Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN (0.01%TFA)
Gradient
Time 0 8 15
B % 5 95 95
Flow Rate: 1.0mL/min
Oven temperature: 40C
Column: SUNFIRE C18, 150*4.6mm, 3.5um

DAD1 A, Sig=214,4 Ref=off



DAD1 B, Sig=254,4 Ref=off



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
9.96	0.04	5760.04	2085.23	99.56
10.18	0.04	25.63	9.88	0.44

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
9.96	0.03	401.29	178.86	100.00