


File D:\DATA\CX4945-03
Acq. Method : D: 1-POS-3MIN.M
Analysis Method : D: 1-POS-3MIN.M
Sample Info : Easy-Access Method: '1-POS-3MIN'
Method Info : Column:SunFire C18 (4.6x 50 mm,3.5 um)
Mobile phase: H2O(0.05%TFA) (A) / ACN(0.05%TFA) (B)
Elution program: Gradient from 5 to100% of B in 1.3 min at 2ml/min
Temperature: 50°C
Detection: UV (214 , 4 nm) and MS (ESI, Pos mode ,110 to 1000 amu)




Integration Results for DAD1 A, Sig=2144 Ref=off

RetTim	Width	Area	Height	Area%
1.38	0.02	5.96	5.02	0.50
1.61	0.01	2232.32	2462.93	99.23
1.81	0.02	5.09	5.13	0.19
2.28	0.02	2.71	2.16	0.08

Integration Results for DAD1 B, Sig=2544 Ref=off

RetTim	Width	Area	Height	Area%
1.61	0.01	1930.88	2252.79	100.00

Ret. Time: 1.61 <<<< POSITIVE SPECTRA >>>>



Ret. Time: 1.81

