



File D:\DATA\Carbenoxolone Sodium-02
Acq. Method : D HPLC
Analysis Method : D HPLC
Sample Info : Easy-Access Method:
Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN (0.01%TFA)
Gradient
Time 0 8 15
B % 5 95 95
Flow Rate: 1.0mL/min
Oven temperature: 40C
Column: SUNFIRE C18, 150*4.6mm, 3.5um



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
10.03	0.04	635.34	246.87	98.52
10.25	0.04	16.15	6.97	1.48

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
8.22	0.03	18.74	8.57	0.70
10.03	0.04	2653.60	1024.88	98.64
10.35	0.04	17.76	6.39	0.66