


Injection Date : 26 Aug 20
 Sample Name : DB07268 -01
 Method Info : Mobile Phase: A:Water(10mM NH₄HCO₃) B:ACN
 Gradient: 5% increase to 95%B within 1.4min, 95%B for 1.6min.
 Flow Rate: 1.8ml/min
 Column: XBridge C18, 4.6*50mm, 3.5um
 Column Temperature: 45 C
 Detection: UV (214, 4nm) and MS (ESI, Neg mode, 110 to 1000 amu)



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
1.51	0.02	1339.49	929.59	57.01
1.51	0.02	1010.17	926.39	42.99

Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
1.52	0.02	2349.99	1575.32	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off


RetTim	Width	Area	Height	Area%
1.52	0.02	1307.53	875.88	100.00

Ret. Time: 1.14

<<<< POSITIVE SPECTRA >>>>



Ret. Time: 1.52



Ret. Time: 1.82

