



File D:\DATA\Deslorelin-LCMS

I
Acq. Method : C:\METHODS\1-POS-3MIN.M
Analysis Method : C:\METHODS\1-POS-3MIN.M
Sample Info : *** No target masses specified ***
Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)
Gradient: 5%-95% B in 1.4min
Flow Rate: 2.3 ml/min
Column: SunFire C18,4.6*50mm,3.5um
Oven Temperature: 50 C
Mass Range: 110-800


*ADC1 A, ELSD




*DAD1 A, Sig=214,4 Ref=off



*DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70



Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
1.34	0.02	492.55	353.41	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off


RetTim	Width	Area	Height	Area%
1.34	0.02	39.83	30.18	100.00

Ret. Time: 0.21


<<<< POSITIVE SPECTRA >>>>




Ret. Time: 0.34




Ret. Time: 0.73




Ret. Time: 1.16




Ret. Time: 1.23




Ret. Time: 1.28




Ret. Time: 1.34



Ret. Time: 1.39



Ret. Time: 1.41



Ret. Time: 1.43

