



File D:\DATA\
Injection Date : 25 Nov 16 10:12 am +0800 Tgt Mass (EZ) :
Sample Name : GDC-0879-04 Location : P2-D-09
Acq. Operator : LCMS Inj : 1
Spec. Reported : MS Integration Inj Volume : 2 ul
Acq. Method : C:\Chem32\1\METHODS\1-POS-3MIN.M
Analysis Method : C:\Chem32\1\METHODS\1-POS-3MIN.M
Sample Info : Easy-Access Method: '1-POS-3min'
Method Info : Mobile Phase: A: water (0.01%TFA) B: ACN (0.01%TFA)
Gradient: 5%-95% B in 1.7 min
Flow Rate: 2.2 ml/min
Column: SunFire C18, 4.6*50mm, 3.5um
Column Temperature: 50 C
MASS Range: 110-900


*ADC1 A, ELSD




*DAD1 A, Sig=214,4 Ref=off



*DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.21	0.02	2207.37	2150.42	99.36
1.27	0.02	17.63	15.68	0.64

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.21	0.01	1647.57	1658.74	100.00

Ret. Time: 1.21 <<<< POSITIVE SPECTRA >>>>

