



Current Data Parameters  
 NAME GSK1120212  
 EXPNO 1  
 PROCNO 1

F2 - Acquisition Parameters  
 Date\_ 20140531  
 Time 12.42  
 INSTRUM spect  
 PROBHD 5 mm PABBO BB-  
 PULPROG zg30  
 TD 65536  
 SOLVENT DMSO  
 NS 16  
 DS 2  
 SWH 8278.146 Hz  
 FIDRES 0.126314 Hz  
 AQ 3.9584243 sec  
 RG 456.1  
 DW 60.400 usec  
 DE 8.00 usec  
 TE 297.6 K  
 D1 1.00000000 sec  
 TD0 1

==== CHANNEL f1 =====  
 NUC1 1H  
 P1 10.70 usec  
 PL1 -3.00 dB  
 SF01 400.1324710 MHz

F2 - Processing parameters  
 SI 32768  
 SF 400.130000 MHz  
 WDW EM  
 SSB 0  
 LB 0.30 Hz  
 GB 0  
 PC 1.00