



Current Data Parameters
NAME Ganetespiib-01
EXPNO 2
PROCNO 1

F2 - Acquisition Parameters
Date_ 20170601
Time 10.29
INSTRUM spect
PROBHD 5 mm PABBO BB-
PULPROG zg30
TD 65536
SOLVENT DMSO
NS 16
DS 2
SWH 8278.146 Hz
FIDRES 0.126314 Hz
AQ 3.9584243 sec
RG 322.5
DW 60.400 usec
DE 8.00 usec
TE 296.8 K
D1 1.00000000 sec
TD0 1

===== CHANNEL f1 =====
NUC1 1H
P1 10.70 usec
PL1 -3.00 dB
SF01 400.1324710 MHz

F2 - Processing parameters
SI 32768
SF 400.1300000 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00

