




Injection Date : 3 Feb 21
Sample Name : Lanifibranor-01
Method Info : Column:SunFire C18 (4.6x 50 mm,3.5 um)
Mobile phase: H2O(0.01%TFA) (A) / ACN(0.01%TFA)(B)
Elution program: Gradient from 10 to95% of B in 1.4 min at 2ml/min
Temperature: 50 °C
Detection: UV (214 , 4 nm) and MS (ESI, Pos mode ,110 to 1000 amu)



Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
1.93	0.03	3769.27	2209.66	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off


RetTim	Width	Area	Height	Area%
1.93	0.03	1635.58	985.37	100.00

Ret. Time: 1.56

<<<< POSITIVE SPECTRA >>>>



Ret. Time: 1.77



Ret. Time: 1.92

