



Agilent Technologies

Sample name: M10126-H-151-03-1

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1D,Sig=230.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
6.965	BB	0.16	7.62	1.78	0.07
7.612	MB m	0.11	3.74	1.16	0.03
7.756	MM m	0.12	3.08	0.85	0.03
8.012	MB m	0.12	6.23	1.92	0.06
9.037	BB	0.12	3.48	1.03	0.03
9.233	BM m	0.36	11246.88	2761.24	99.28
11.540	MM m	0.17	25.80	5.46	0.23
11.633	MM m	0.14	16.16	3.75	0.14
12.602	MM m	0.16	15.29	3.42	0.13