



Agilent Technologies

Sample name: M10152-Naporafenib-02

Instrument: A012

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1A, Sig=220.0, 4.0 Ref=off

| RT [min] | Peak | Area   | Height | Area%  |
|----------|------|--------|--------|--------|
| 5.745    |      | 10.8   | 3.6    | 0.116  |
| 6.885    |      | 10.2   | 2.8    | 0.109  |
| 8.215    |      | 7.7    | 2.3    | 0.082  |
| 8.346    |      | 3.2    | 1.0    | 0.034  |
| 8.584    |      | 4.4    | 1.0    | 0.047  |
| 9.041    |      | 9314.1 | 1969.0 | 99.612 |