



Agilent Technologies

Sample name: M10192-Seladelpar sodium salt-01

Injection date:

Instrument: 5

Location:

Inj. volume: HPLC


Type:

Acq. method:

Sample amount:

Processing method:

Manually modified:



Signal: DAD1D,Sig=230,4 Ref=off

RT [min]	Type	Width [min]	Area	Height	Area%
8.230	MB m	0.12	6.51	1.96	0.03
8.633	MM m	0.09	2.25	0.75	0.01
8.973	MB m	0.13	9.94	2.68	0.04
9.104	BB	0.11	7.14	2.30	0.03
10.987	MB m	0.14	5.96	1.57	0.03
11.825	BB	0.84	21985.15	2695.14	99.32
12.068	BB	0.16	25.53	6.05	0.12
12.211	BB	0.14	11.60	3.14	0.05
12.403	BB	0.23	13.06	2.26	0.06
13.174	MB m	0.12	3.88	0.94	0.02
13.458	MB m	0.16	12.91	2.76	0.06
13.909	MB m	0.21	42.13	8.80	0.19
14.788	MB m	0.14	8.50	2.18	0.04