



Sample name: M10232-Cefotaxime-02

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1F,Sig=260.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
2.284	BB	0.19	15.03	5.37	0.20
2.537	MM m	0.15	1.67	0.46	0.02
4.210	MM m	0.11	1.28	0.36	0.02
4.312	MB m	0.15	20.17	6.05	0.27
4.513	MM m	0.11	0.91	0.21	0.01
4.731	MM m	0.21	14.27	3.29	0.19
4.829	MM m	0.13	1.24	0.34	0.02
5.077	BB	0.21	22.88	7.27	0.30
5.258	BV	0.32	7435.68	2300.31	98.61
5.466	VM m	0.06	7.62	2.33	0.10
5.520	MM m	0.07	5.13	2.11	0.07
5.586	MB m	0.14	2.13	0.48	0.03
5.788	BB	0.13	4.39	1.46	0.06
6.226	BM m	0.20	6.40	1.66	0.08
7.199	MM m	0.15	1.67	0.52	0.02