

Analysis Report

<Sample Information>


Sample Name : M10671-Garcinone E-01
Instrument : SHIMADZU LC-20AT
Column : Wonda Cract ODS-2 (5 µm, 4.6×250 mm)
Column Temp. : 35°C
Mobile phase flow : 1.0 ml/min
Injection Volume : 5 µL
Detection Wave. : UV-220nm
Con. & Solvent : 0.2mg/ml (Methanol)

<Analysis condition>

50%→100% Acetonitrile in 0.05% Phosphoric acid H2O 13 min
100%→100% Acetonitrile in 0.05% Phosphoric acid H2O 3 min

<Chromatogram>

mAU



<Peak Table>

PDA Ch1 220nm

Peak#	Ret. Time	Area	Height	Conc.	Unit	Area%
1	12.539	13701	1747	0.000		0.373
2	12.903	3619873	731216	0.000		98.496
3	13.048	41557	11978	0.000		1.131
		3675131	744942			100.000