



Agilent Technologies

Sample name: M10871-PF-07104091-02

Injection date:

Instrument: A010

Location:

Inj. volume: 3


Type:

Acq. method: HPLC

Sample amount:

Processing method:

Manually modified:



Signal: DAD1F, Sig=260.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
6.082	BM m	0.19	9.51	2.53	0.07
6.356	MM m	0.18	0.65	0.14	0.00
6.630	MM m	0.10	0.51	0.17	0.00
6.753	BM m	0.37	14140.22	2980.29	99.61
7.054	MM m	0.07	10.14	2.83	0.07
7.099	MM m	0.16	8.91	2.09	0.06
7.465	BB	0.17	10.45	2.47	0.07
7.647	MM m	0.19	2.32	0.41	0.02
9.145	MM m	0.72	13.26	0.61	0.09