



Agilent Technologies

Sample name: M10927-GFB-8438-01

Instrument: A010

Inj. volume: 2

Acq. method: HPLC

Processing method:


Manually modified:

Injection date:

Location:

Type:

Sample amount:



Signal: DAD1D,Sig=230.0,4.0 Ref=360.0,100.0

RT [min]	Type	Width [min]	Area	Height	Area%
6.832	MM m	0.12	10.04	3.00	0.13
7.175	MM m	0.08	1.67	0.60	0.02
7.335	BM m	0.38	7901.72	2233.02	99.61
7.772	MB m	0.12	11.47	3.74	0.14
8.579	MM m	0.11	1.79	0.55	0.02
9.033	MM m	0.08	2.06	0.74	0.03
9.389	MB m	0.14	4.28	1.02	0.05